

# NobelReplace<sup>®</sup> Conical Connection

## Product overview

### High initial stability, even in compromised bone situations

The tapered implant body is designed for high initial stability and is ideal for Immediate Function in both extraction sockets and healed sites.

### Strong sealed connection

Advanced internal conical connection with hexagonal interlocking offers a tight seal and high mechanical strength.



### Broad range of restorative solutions

Prefabricated and CAD/CAM NobelProcera restorations to support all temporary and final solutions.

### Natural-looking esthetics


Built-in platform shifting designed to increase soft tissue volume.

### Enhanced osseointegration

Unique oxidized TiUnite surface with grooves increases implant stability through faster bone formation and ensures long-term success.




### Efficient treatment flow

Step-by-step drilling protocol simplifies site preparation and ensures predictable outcomes in all indications.



# Order your implant of choice today.




## NobelReplace Conical Connection

Platform	Implant Ø	Length				
		8mm	10mm	11.5mm	13mm	16mm
	3.5mm	36699	36700	36701	36702	36703
	4.3mm	36704	36705	36707	36708	36709
	5.0mm	36710	36711	36712	36713	36714



With TiUnite on collar

## NobelReplace Conical Connection PMC (Partially Machined Collar)\*


Platform	Implant Ø	Length				
		8mm	10mm	11.5mm	13mm	16mm
	3.5mm	37284	37285	37287	37288	37289
	4.3mm	37290	37291	37292	37293	37294
	5.0mm	37295	37296	37297	37298	37299



With 0.75 mm machined collar

\* Cover screw included.

## Cover screw

	NP	RP
	36649	36650




On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit [nobelbiocare.com/warranty](http://nobelbiocare.com/warranty)



For the full assortment of Nobel Biocare products, see our current product catalog.

[nobelbiocare.com/catalog](http://nobelbiocare.com/catalog)



### Order online

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

[nobelbiocare.com/store](http://nobelbiocare.com/store)



### Order by phone or fax

Call our customer service, send a fax or contact your sales representative.



# Prosthetic components for implants with conical connection.

## Temporary abutments<sup>1</sup>

		NP			RP		
	<b>Healing Abutment</b>	<b>3 mm</b>	<b>5 mm</b>	<b>7 mm</b>	<b>3 mm</b>	<b>5 mm</b>	<b>7 mm</b>
	Ø3.6	36639	36640	36867	36643	36644	36872
	Ø5.0	36641	36642	36868	36645	36646	36873
	Ø6.0	-	-	-	36647	36648	36874
	<b>Immediate Temporary Abutment<sup>2</sup></b>	<b>1.5 mm</b>	<b>3 mm</b>		<b>1.5 mm</b>	<b>3 mm</b>	
		36653	36655		36654	36656	
	<b>Temporary Abutment Engaging</b>	36663			36664		
	<b>Temporary Coping Multi-unit Titanium<sup>4</sup></b>	29046			29046		

		NP			RP		
	<b>Healing Abutment Bridge</b>	<b>3 mm</b>	<b>5 mm</b>	<b>7 mm</b>	<b>3 mm</b>	<b>5 mm</b>	<b>7 mm</b>
	Ø4.0	36864	36865	36866	-	-	-
	Ø5.0	-	-	-	36869	36870	36871
	<b>QuickTemp<sup>TM</sup> Abutment<sup>3</sup></b>	<b>1.5 mm</b>	<b>3 mm</b>		<b>1.5 mm</b>	<b>3 mm</b>	
		36659	36657		36660	36658	
	<b>Temporary Abutment Non-Engaging</b>	36661			36662		
	<b>Temporary Coping Multi-unit Plastic<sup>5</sup></b>	DCA			DCA		
		468-0			468-0		

## Impression copings

### Implant level

		NP 13 mm	RP 13 mm	RP (low profile) 9 mm
	<b>Closed Tray</b>			
	Ø3.6	36538	36540	36541
	Ø5.0	36539	36542	36543
	Ø6.0	-	36544	36545

		NP 10 mm	NP 14 mm	RP 10 mm	RP 14 mm
	<b>Open Tray</b>				
	Ø3.6	36258	36260	36263	36262
	Ø5.0	36259	36261	36265	36264
	Ø6.0	-	-	36267	36266

		NP	RP
	<b>Bridge Open Tray</b>	36930	36931


### Abutment level

		NP	RP
	<b>Closed Tray Multi-unit</b>	29090	29090
	<b>Bar Closed Tray Multi-unit</b>	29093	29093
	<b>Snappy<sup>TM</sup> Abutment</b>		
	4.0	35995	35996 35997 Wide
	5.5	35868	35869 35870 Wide

		NP	RP
	<b>Open Tray Multi-unit<sup>6</sup></b>	29089	29089
		<b>4/pkg</b>	<b>20/pkg</b>
	<b>Locator<sup>®</sup> Abutment</b>	REF08505	REF08505-20
	<b>Block Out Spacer</b>	-	REF08514
		<b>2/pkg</b>	<b>10/pkg</b>
	<b>Locator<sup>®</sup> Male Processing Package</b>	REF08519-2	REF08519-10

## Prefabricated final abutments<sup>1</sup>

### Zirconia

	NP	RP
	<b>Procera® Esthetic Abutment</b>	
	NP #1 36627	RP #5 36631
	NP #2 36628	RP #6 36632
	NP #3 36629	RP #7 36633
	NP #4 36630	RP #8 36634
		RP #9 36635
		RP #10 36636
		RP #11 36637
		RP #12 36638

**Consistent color-coding**  
Implants, surgical and prosthetic components are color-coded for accurate and fast identification.









Narrow Platform NP




Regular Platform RP


### Titanium

	NP		RP	
	<b>Esthetic Abutment</b>			
	1.5 mm 3mm	36665 36666	1.5 mm 3mm	36669 36671
	4.5 mm	36249	4.5 mm	36251
	<b>Snappy™ Abutment 4.0°</b>			
	1.5 mm 3mm	36695 36696	1.5 mm 3mm	36693 36694
				36691 Wide 36692 Wide
	<b>Narrow Profile Abutment</b>			
	7 mm 9 mm	36678 36679	7 mm 9 mm	36680 36681



	NP		RP	
	<b>Esthetic Abutment 15°</b>			
	1.5 mm 3mm	36667 36668	1.5 mm 3mm	36672 36673
	4.5 mm	36250	4.5 mm	36252
	<b>Snappy™ Abutment 5.5°</b>			
	1.5 mm 3mm	36684 36688	1.5 mm 3mm	36682 36683
				36689 Wide 36690 Wide
	<b>Locator® Abutment</b>			
	1.0 mm	REF02069		REF02075
	2.0 mm	REF02070		REF02076
	3.0 mm	REF02071		REF02077
	4.0 mm	REF02072		REF02078
	5.0 mm	REF02073		REF02079
	6.0 mm	REF02074		REF02080





### Gold

	NP	RP
	<b>GoldAdapt™ Engaging</b>	
	36728	36729

	NP	RP
	<b>GoldAdapt™ Non-Engaging</b>	
	36726	36727

### Multi-unit

	NP	RP
	<b>Multi-unit Abutment</b>	
	1.5 mm	36611 36616
	2.5 mm	36613 36617
	3.5 mm	36624 36625
	4.5 mm	- 36626
	<b>Healing Cap Multi-unit</b>	
	1/pkg	31145 31145
	5/pkg	29064 29064
	Wide 1/pkg	31146 31146

	NP		RP	
	<b>17° Multi-unit Abutment</b>			
	2.5 mm 3.5 mm	36614 36615	2.5 mm 3.5 mm	36618 36619
	<b>30° Multi-unit Abutment</b>			
	3.5 mm 4.5 mm	36620 36621	3.5 mm 4.5 mm	36622 36623
	<b>Gold Coping Multi-unit<sup>4</sup></b>			
	1/pkg	29043		29043
	5/pkg	29042		29042
	<b>Gold Coping Bar Multi-unit<sup>4</sup></b>			
	29045			29045

1 Abutment screw included with all abutments.

2 Plastic coping 31656 included (cannot be used as a burn-out).

3 Plastic coping 33404 included (cannot be used as a burn-out).

4 Prosthetic screw included.

5 Prosthetic screw not included (cannot be used as a burn-out).

6 15 mm guide pin included.

7 Includes abutment screw, impression coping, healing cap and plastic/temporary coping engaging (plastic/temporary coping non-engaging also available).



### NobelProcera® individualized CAD/CAM restorations



- Individualized abutments, screw-retained crowns and implant bridges in zirconia (4 shades) and titanium
- Fixed and fixed-removable implant bars overdenture
- Coping and bridge frameworks and full-contour restorations in a variety of materials



Contact your local Nobel Biocare sales office to learn more about the full range of precision-milled NobelProcera restorations.






For all Nobel Biocare and other major implant systems.




### Replicas

	NP	RP		
 <b>Implant</b>	36697	36698		
 <b>Snappy™ Abutment</b>				
<b>4.0</b>	35982	35983	35984 Wide	
<b>5.5</b>	35854	35855	35856 Wide	




	NP	RP		
 <b>Multi-unit Abutment</b>	<b>1/pkg</b> 31161	<b>5/pkg</b> 29110	<b>1/pkg</b> 31161	<b>5/pkg</b> 29110
 <b>Locator® Abutment</b>		<b>4/pkg</b>	<b>20/pkg</b>	
	<b>Ø4 mm</b>	REF08530	REF08530-20	
	<b>Ø5 mm</b>	REF08516	REF08516-20	



### Screws

	NP	RP		
 <b>Abutment Screw</b> (for zirconia and titanium restorations)	36917	36918		
 <b>Abutment Screw Multi-unit Angled</b>	36892	36893		
 <b>Prosthetic Screw Multi-unit</b>	29285	29285		

	NP	RP		
 <b>Laboratory Screw</b> only for titanium restorations	31168	29290 (5/pkg)		
 <b>Laboratory Screw Multi-unit Angled</b>	31164	31165		
 <b>Laboratory Screw Multi-unit</b> 5/pkg	29287	29287		

### Other components

	NP	RP		
 <b>Protection Analog</b> 5/pkg	36730	36731		
 <b>32309 Prosthetic Kit</b>				
Screwdriver Machine Unigrip™ 20mm	29151			
Screwdriver Machine Unigrip™ 30mm	29153			
Screwdriver Machine Multi-unit 21mm	29158			
Manual Torque Wrench Prosthetic	29165			

	NP	RP		
 <b>Protection Analog Multi-unit</b> 5/pkg	29123	29123		
 <b>Guide Pin Multi-unit</b>				
10mm (1/pkg)		31154		
10mm (5/pkg)		29102		
20mm (1/pkg)		31155		
20mm (5/pkg)		29103		

# Surgical components.

Precision Drill	
NP, RP	36118



Drill with Tip Tapered Ø 2 mm	
NP, RP	36117



Drills	8 mm	10 mm	11.5 mm	13 mm	16 mm
NP 3.5 mm	32075	29367	36113	29368	29369
RP 4.3 mm	32076	29370	36114	29371	29372
WP 5.0 mm	32077	29373	36115	29374	29375

Dense Bone Drills	13 mm	16 mm
NP 3.5 mm	29377	29378
RP 4.3 mm	29380	29381
WP 5.0 mm	29383	29384



Screw Taps	
NP 3.5 mm	36717
RP 4.3 mm	32090
WP 5.0 mm	32091



Implant Drivers	28 mm	37 mm
NP	36718	36719
RP	36720	36721



Bone Mills and Guides	Bone Mill with Guide	Bone Mill Guide
NP	36722	36723
RP	36724	36725



**36894 NobelReplace Tapered Surgery Kit**  
Includes instruments and drills to perform surgery for all NobelReplace/Replace Select Tapered implants Ø 3.5, 4.3 and 5.0 mm.



**36915 Implant Driver Kit Conical Connection**  
Additionally needed for inserting NobelReplace Conical Connection implants.



## The key to successful treatments

Diagnose and plan implants in 3D based on the patient's anatomy and prosthetic needs with the NobelClinician Software. For more complex cases, you have the option to place the implants fully guided with the well-established NobelGuide concept.

[nobelbiocare.com/nobelclinician](http://nobelbiocare.com/nobelclinician)



[nobelbiocare.com/nobelreplacecc](http://nobelbiocare.com/nobelreplacecc)